
SPECTROPHOTOMETERS


S-200 S-220



MAIN MENUE



LARGE DISPLAY MEASURE



TIME SCAN RESULTS AND DATA PROCESSING

BOECO SPECTROPHOTOMETER MODELS S-200 VIS & S-220 UV/VIS

The BOECO S-220 (UV/VIS) and S-200 (VIS) are high quality, compact, low cost measurement systems for daily analysis in education, QC and basic research.

▶ Compact single beam optics with full range scanning

The single beam optics are compact and bench space saving. The long life Hamamatsu Xenon lamp optics in the S-220 ensure quick and reliable performance and the Tungsten Halogen lamp used in S-200 also provide a reliable measurement.

▶ Color touch screen operation

The intuitive color touch screen operation provides simple access to an extensive range of functions. The touch screen is sensitive to stylus and laboratory gloves. Icon driven on-board software improves accessibility and the graphical display allows spectrum or standard curve to be shown on the screen. The forward and back quick key allows the user to proceed or swiftly return to the process. An enlarged data display for photometry measurement makes result reading easier.

▶ Various measurement modes

Operation modes include photometric, multiple wavelength analysis, spectrum scanning, time scan and kinetics; direct concentration results are included.

Optional accessories

A variety of accessories are included such as test tube holder, flow cell with sipper, temperature control holder, long path length cuvette holder & multiple cell holder are available to enhance different application needs.

▶ Storage and data output

External storage with SD card and free downloadable PC Software MasterReport (www.boeco.com) allows data export to PC in compatible text or spreadsheet format for further data processing in the PC. Method and result storage is almost unlimited by exchanging SD card when needed. Printer options are available for direct result printing with graphics.

▶ Validation function

To ensure optimum instrument performance, self diagnosis functions are equipped in GLP/GMP feature for performance validation and auditing

Code	Description
BOE 8620000	Model S-200 Vis Spectrophotometer, single beam
	with full range scanning and color touch screen
	operation. Supplied with 10 x 10 mm cuvette holder
BOE 8622000	Model S-220 UV/Vis Spectrophotometer, single beam
	with full range scanning and color touch screen
	operation. Supplied with 10 x 10 mm cuvette holder

Specification

Wavelength Range:

Spectral Bandwidth:

Transmittance accuracy: Transmittance repeatability:

Baseline flatness: Noise level: Baseline stability:

Stray light:

Wavelength controlled variable: Wavelength accuracy: Wavelength repeatability: Wavelength scan speed:

Wavelength move speed:

Absorbance: Transmittance: Spectrum Scanning: Concentration: Selectable Resolution:

Light source: Detector:

Display screen: Printer:

Metering mode: Memory: Time Scan: Analysis: GLP: Size:

Power requirement: Power consumption: Communication ports:

Weight:

Code

S-200 Vis

320 to 1100nm

6nm

±0,5% T (NIST 930 Filter)

0.2% T

± 0.002 Abs (330-1090nm) $\leq 0.001 \text{ Abs } (500 \text{nm})$ $\leq 0.001 \text{ Abs/h } (500 \text{nm})$

(after 2 hours warm up)

≤ 0.5% T 0.2nm + 1nm ≤ 0.5nm 2400nm/min

(0,2 sampling interval without filter) to any specified position within 1sec.

-0.3 to 1.999 0 to 199.9%

Yes

-300 to 1999

1, 0.1, 0.01 or 0.001 Tungsten Halogen lamp S-220 UV/Vis

190 to 1000nm

5nm

±1% T (NIST 930 Filter)

0.5% T

± 0.005 Abs (200-990nm) $\leq 0.005 \text{ Abs } (250 \text{ nm})$

 $\leq 0.005 \text{ Abs/h } (250 \text{nm})$

(after 2 hours warm up) ≤ 0.5% T

0,2nm + 2nm ≤ 1nm 300nm/min

(0,2 sampling interval without filter) to any specified position within 1sec.

-0.3 to 1.999 0 to 199.9% Yes

-300 to 1999 1, 0.1, 0.01 or 0.001

pulsed-Xenon lamp

Silicon photodiode 4,3 inches colorful touch LCD screen specified 80-column thermal printer (series port)

> Single beam SD card storage

Graphical and calculated concentration value Absorbance and wavelength of peaks and valleys Real time clock and calendar, Self Diagnosis 400 (W) x 280 (D) x 160 (H) mm AC, 100-240V, 50/60Hz

100VA

Serial printer port connects thermal printer USB port connects PC

SD card port saves data and measurement methods Accessories port connects and controls serval options 4 kgs

Accessories

Test tube holder (only for S-200)
Rectangular long-path cuvette holder for cuvettes
with 10, 20, 30, 50 and 100 mm path-length
Micro-cuvette holder, for cuvettes with centre
height of 15 mm
Flow cuvette holder, incl. quartz glass flow cuvette
of 150 µl
Set of Auto sample sipper and Flow cuvette holder
with quartz glass flow cuvette of 150 μl
Electronic thermostat (Peltier element) TC cuvette
holder (only for S-220, S-300)
Automatic 5 position cuvette holder
Thermo printer with 100V-240V AC power supply
Tungsten halogen lamp (S-200)
Xenon Lamp module (S-220, S-300)
UV DETECTIVE software to control and operate the spectrophotometer on a PC. The versatile software can control all spectrophotometer operations such as photometry, wavelength scans, time scans and more Further functions include storage of methods programs, saving of numerical and graphical data, downstream data processing, data transfer to commercial spreadsheets such as Excel® and report generation

Description



TEST TUBE HOLDER



5-PLACE AUTO SAMPLE HOLDER



FLOW CUVETTE HOLDER



THERMO PRINTER



LONG PATH **CUVETTE HOLDER**



ELECTRONIC THER-MOSTAT HOLDER



SAMPLE SIPPER



MICRO CUVETTE HOLDER



S-300



MAIN MENUE





BOECO LIFE SCIENCE SPECTROPHOTOMETER MODEL S-300

The BOECO S-300 life science spectrophotometer allows measurement of nucleic acid concentrations and purity (using ratio function) including protein concentrations. As a high quality spectrophotometer, the S-300 features touch screen operation packaged as a lightweight system with a compact footprint for life science and education related applications.

▶ Life Science Programs

The S-300 contains onboard functions for the quantification of nucleic acid, including dsDNA, ssDNA, RNA and Oligonucleotides. The purity of the nucleic acid can also be determined with the ratio A260/A280 calculation. Protein concentrations can be measured from a range of colourimetric assays such as Bradford, Lowry, Biuret and BCA. Standard calibration data and curves can also be displayed. Furthermore, proteins can be quantified at 280nm. Bacterial cell density at 600nm can also be measured under the OD600 cell culture optical density function. It can define a bacterial culture in exponential growth phase and at the most appropriate time for harvest or induction.

▶ Compact Optics with Full Range Scanning

The single beam optics are compact resulting in significant bench space saving. The long life Hamamatsu Xenon lamp optics system in the S-300 ensures quick and reliable performance.

▶ Color Touch Screen Operation

The intuitive color touch screen provides simple access to an extensive range of function. The touch screen is sensitive to stylus or hands (with and without gloves). Icon driven on board software improves accessibility and the quick action keys are another convenience feature.

Various measurement modes

In addition to the Lifescience program, the S-300 also features conventional spectrophotometer functions such as single/ multiple wavelength analysis, spectrum scanning, kinetics and concentration measurement.

Optional accessories

A various selection of optional accessories is available such as flow cell with sipper, temperature control holder, long path length cuvette holder & multiple cell holder to enhance different application needs.

▶ Storage and data output

External storage with SD card allows data export to PC in compatible text or spreadsheet format. Free downloadable PC Software MasterReport (www. boeco.com) allows data export to PC in compatible text or spreadsheet format for further data processing in the PC.

Method and result storage is almost unlimited by exchanging SD card when needed. Printer options are available for direct result printing with graphics.

▶ Validation function

To ensure optimum instrument performance, self diagnosis functions are equipped in GLP/GMP feature for performance validation and auditing.

Code

Description

BOE 8630000

Model S-300 UV/Vis Life Science Spectrophotometer, single beam with full range scanning and color touch screen operation. Supplied with installed micro cuvette holder (centre height 15 mm) and with optional 10 x 10 mm cuvette holder and sample pack of 8 pcs. disposable UV Micro Cuvettes

Specification S-300 UV/Vis

Wavelength Range: 190 to 1000nm

Wavelength Resolution 0.2nm Spectral Bandwidth: 5nm

Transmittance accuracy: $\pm 1\%$ T (NIST 930 Filter)

Transmittance repeatability: 0.50% T

Detection limit Concentration: dsDNA 1.5 - 100µh/ml (for 100µl cell)

Noise level: 0.005 Abs (at 250 nm) Stray light: ≤ 0.5% T at 220, 340 nm

Wavelength accuracy: ± 2nm

Wavelength repeatability: ≤ 1nm

Absorbance: -0.3 to 1.999

Transmittance: 0 to 199.9%

Spectrum Scanning: Yes

Concentration: 0 to 1999

Light source: pulsed-Xenon lamp Detector: Silicon photodiode

Display screen: 4,3 inches colorful touch LCD screen

Printer: specified 80-column thermal printer (series port)

Metering mode: Single beam Memory: SD card storage

Time Scan: Graphical and calculated reaction activity

Wavelegth Scan Analysis: Absorbance and wavelength of peaks and valleys GLP: Real time clock and calendar, Self Diagnosis

Size: 400 (W) x 280 (D) x160 (H) mm Power requirement: AC, 100-240V, 50/60Hz

Power consumption: 700VA

Communication ports: Serial printer port connects thermal printer

USB port connects PC

SD card port saves data and measurement methods Accessories port connects and controls serval options

Weight: 4 kgs

TRAYCELL

The HELLMA® TrayCell is a fibre-optic ultra-micro cell designed to the UV/Vis analysis of DNA/RNA and proteins. The dimensions of the TrayCell are equivalent to a standard cuvette in order to work in most sprectrophotometers.

▶ Efficient accessory

for your spectrophotometer

▶ Extremely flexible and cost-effective solution

for the analysis of very small sample volumes (0,7 - 5µl)

Ideal for biomolecular laboratories

to perform the analysis of nucleic acids and proteins in very small volumes

Specification TrayCell 105.810-UVS

Window material: Quartz SUPRASIL Wide/depth 12,5 x 12,5 mm Height: 59,5 mm Volume: 0,7 - 5 μ l

Light path: 0,2 mm or 1 mm (\pm 0,02) depending on the cap

Max. temperature: 50 °C Centre height: 15 mm

Fibre optic cable: built-in, not exchangeable

UV/Vis low solarisation 190 nm - 1.100 nm

Code Description

HEL 105810-A3-V1-46	TrayCell 105.810-UVS, centre height 15 mm
HEL 665-703-1-40	TrayCell Cap 665.703, 1 mm Light path
HEL 665-704-0.2-40	TrayCell Cap 665.704, 0,2 mm Light path













HEL 100-OS/10



HEL 105.202-QS

CUVETTES

Disposable Cuvettes:

Code	Description
BRA 759007	Disposable Macro Cuvettes, PS,
	Window: 10 x 35 mm; 10 mm light path
	Filling volume: min 2,5 / max 4,5 ml
	grouped by mold cavity number, neutral packing
	Wavelength: From 340 to 900 nm
	Packing: 100 pcs./box, 1000 pcs./carton
BRA 759017	Disposable Semi-Micro Cuvettes, PS,
	Window: 4,5 x 23 mm; 10 mm light path
	Filling volume: min 1,5 / max 3,0 ml
	grouped by mold cavity number, neutral packing
	Wavelength: From 340 to 900 nm
	Packing: 100 pcs./box, 1000 pcs./carton
BRA 759170	BRAND Disposable UV Macro Cuvettes,
	Filling volume: min 2,5 / max 4,5 ml
	grouped by mold cavity number,
	Wavelength: From 220 to 900 nm,
	10 mm light path, pack of 100 pcs.
BRA 759150	BRAND Disposable UV Semi-Micro Cuvettes,
	Filling volume: min 1,5 / max 3,0 ml
	grouped by mold cavity number,
	Wavelength: From 220 to 900 nm,
	10 mm light path, pack of 100 pcs.
BRA 759220	BRAND Disposable UV Micro Cuvettes,
	Center height 15 mm, Vol. 70 µl up to 550 µl,
	Wavelength: From 220 to 900 nm,
	10 mm light path, pack of 100 pcs.

Optical and quartz glass Cuvettes: