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DS 450 Demineraliser

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Operating instructions

Bedienungsanleitung für DS450, Art. Nr.: 11042506/07.

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General notes and safety precautions

- Read these notes and precautions for your own well-being and proper functioning of the cartridge!
- The deionized water that is produced is not fit for drinking.
- Use only the original raw water hose, as the flow restrictor integrated in it prevents damage to the cartridge.
- **Never** insert a shut-off device in the cartridge outlet nor kink the outlet hose, as the plastic cartridge could burst under the resulting pressure build up. Ensure free water flow at all times.
- The use of this system in any other way than that described in these operating instructions invalidates the guarantee.
- An earthing contact socket outlet is required for the electrical connection.
- Under some circumstances, improper handling could result in the leakage of ion exchange resin from an ion exchange cartridge. **CAUTION: Danger of slipping** on leaked out resin! As the resin material can cause slight irritation on contact with skin and/or eyes, immediately turn off the raw water tap to stop further leakage and **proceed as follows**: on skin contact, immediately wash contacted skin areas with plenty of water. On contact with eyes, remove contact lenses if this is possible without difficulty, then immediately and carefully wash eyes with plenty of water for about 10 minutes. When irritation still occurs, seek medical advice immediately.
- Wear protective gloves to collect up the leaked out resin, fill it into a clearly labelled plastic container and return it to the authorized specialist supplier for waste disposal.

Mounting and assembling instructions

1. Draw the cartridge out from the bracket. Mount it near to the raw water tap using the enclosed screws and dowels. Ensure that the DS 450 mixed-bed ion exchanger is freely accessible

CAUTION: Before starting step 1, please make sure that there are no electrical wires or water-carrying pipes in the vicinity of the planned drilling holes.

2. After mounting the bracket, use slight pressure to insert the cartridge with attached conductivity meter in it.
3. Connect the raw water hose, free of kinks and coils, to the water tap and to the R 3/4" thread at the bottom of the cartridge.

CAUTION: The straight connector of the raw water hose contains the water restrictor and must be connected to the water tap. The angled connector is to be fitted to the R 3/4" thread of the cartridge.

4. Fit the pure water hose to the hose nipple of the conductivity meter.

CAUTION: Ensure that the pure water hose is neither blocked in any way nor kinked, as the plastic cartridge could then burst under the resulting build-up of pressure.

5. Connect the power plug to the appropriate socket with earth contact: 115V/230V, 50/60 Hz.
6. Carefully open the water tap and watch for any leaks. Should there be leakage at the union nut of the conductivity meter or at the hose connections, turn the water tap off and tighten the union nut or hose connections by hand until the water leak is stopped.
7. The system is now ready to produce demineralized water. Keep a watch on the conductivity meter, as effluent water is first to be taken for use when the pointer of the conductivity meter is in the green section or at a maximum of 20 µS/cm.

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Change of cartridge

1. Close the water tap and remove the electric plug from the socket, then disconnect the raw water hose from the cartridge.

CAUTION: When opening, water may leak out of the cartridge. Make sure, that this does not impair any other systems in the vicinity.

2. Remove the cartridge with conductivity meter from the bracket.
3. Unscrew the conductivity meter from the exhausted cartridge.

4. Fit the conductivity meter on the fresh cartridge and fit the combination in the bracket. Connect the raw water hose to the fresh cartridge. Only hand tighten the hose nipple.

5. Plug in the power plug. The system is now ready for use again. Please note: When no water flows through the cartridge, re-ionization effects could cause the conductivity meter to show incorrect values. Please do not use the water produced until the flow of it through the measuring cell brings the conductivity meter display to a correct value.

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

Applicable EEC Directives:
 EEC- Directive low voltage
 EEC- Directive electromagnetic compatibility
 EMV 2004/108/EG

Applied standard:
 55011 und 61000

Niederahr, June 2014

Manufacturer/ CE-Representative: Leo Trumm

Note on the waste disposal of equipment

According to your state government requirements and the 2002/96 EC and 2006/66/EC directives, equipment that is to be scrapped can be brought to authorized collection points for recycling. Alternatively, it can be returned to us for proper recycling/waste disposal. In the case of return for repair, incorrect delivery or double delivery, please use the original cardboard box/packaging whenever possible. Send fall and knock protected.

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Boeckel + Co (GmbH + Co)