

# JIPO PORCELAIN

## CRUCIBLES, LOW FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	19	5	JIP 100301	5
35	22	10	JIP 100351	5
40	25	17	JIP 100401	5
45	28	21	JIP 100451	5
50	32	34	JIP 100501	5
60	38	62	JIP 100601	5
70	44	91	JIP 100701	5



## CRUCIBLES, MIDDLE FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	25	10	JIP 200301	5
35	28	12	JIP 200351	5
40	32	20	JIP 200401	5
45	36	30	JIP 200451	5
50	40	45	JIP 200501	5
60	48	80	JIP 200601	5
70	56	120	JIP 200701	5
80	82	200	JIP 200801	5



## CRUCIBLES, HIGH FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	38	15	JIP 300301	5
35	44	26	JIP 300351	5
40	50	35	JIP 300401	5
45	56	50	JIP 300451	5
50	62	72	JIP 300501	5
60	75	120	JIP 300601	5



## LIDS FOR CRUCIBLES

Ø (mm) crucible	Ø (mm) lid	Code	Pack of
30	34	JIP 000301	5
35	39	JIP 000351	5
40	44	JIP 000401	5
45	49	JIP 000451	5
50	54	JIP 000501	5
60	64	JIP 000601	5
70	74	JIP 000701	5
80	84	JIP 000801	5



## GOOCH CRUCIBLES

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	36	17	JIP 234101	5
35	43	25	JIP 234301	5
39	49	35	JIP 234401	5
60	71	120	JIP 234501	5



## FILTER CRUCIBLES

with porous base



H (mm)	Vol (ml)	Porosity	Code	Pack of
28	8	P1 = 6 $\mu$	JIP 400256	2
28	8	P2 = 7 $\mu$	JIP 400257	2
28	8	P3 = 8 $\mu$	JIP 400258	2
35	15	P1 = 6 $\mu$	JIP 400306	2
35	15	P2 = 7 $\mu$	JIP 400307	2
35	15	P3 = 8 $\mu$	JIP 400308	2
40	25	P1 = 6 $\mu$	JIP 400356	2
40	25	P2 = 7 $\mu$	JIP 400357	2
40	25	P3 = 8 $\mu$	JIP 400358	2
43	35	P1 = 6 $\mu$	JIP 400406	2
43	35	P2 = 7 $\mu$	JIP 400407	2
43	35	P3 = 8 $\mu$	JIP 400408	2
50	50	P1 = 6 $\mu$	JIP 400456	2
50	50	P2 = 7 $\mu$	JIP 400457	2
50	50	P3 = 8 $\mu$	JIP 400458	2

## BUECHNER FUNNELS



Nom. size	Stem length	Hole $\varnothing$	Code	Pack of
45	43	1	JIP 237101	1
55	64	1	JIP 237201	1
70	64	1	JIP 237301	1
90	71	2	JIP 237401	1
110	83	2	JIP 237501	1
125	85	2	JIP 237601	1
150	96	2	JIP 237701	1
185	117	2	JIP 237801	1
240	117	2	JIP 237901	1
270	130	2	JIP 237921	1
320	133	2	JIP 237991	1

## DESICCATOR PLATES, ACC. TO CSN



$\varnothing$ (mm)	for Desic.	Code	Pack of
90	100	JIP 238101	1
140	150	JIP 238331	1
189	200	JIP 238521	1
241	250	JIP 238821	1
290	300	JIP 238991	1

## DESICCATOR PLATES, ACC. TO DIN



$\varnothing$ (mm)	for Desic.	Code	Pack of
90	100	JIP 138101	1
140	150	JIP 138201	1
190	200	JIP 138301	1
235	250	JIP 138401	1
280	300	JIP 138501	1