

EVAPORATION DISH, ROUND BOTTOM

deep form, with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
40	16	10	JIP 206001	5
50	20	20	JIP 206101	5
63	25	30	JIP 206121	5
80	32	60	JIP 206201	5
81	38	107	JIP 206221	5
94	42	140	JIP 206301	5
100	40	150	JIP 206321	5
112	47	232	JIP 206401	5
125	50	285	JIP 206411	5
132	55	300	JIP 206421	5
138	61	400	JIP 206501	5
160	64	580	JIP 206611	5
170	74	700	JIP 206621	1
200	80	1000	JIP 206721	1
202	86	1300	JIP 206801	1
320	140	6000	JIP 206991	1



EVAPORATION DISH, ROUND BOTTOM

semideep form, with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
72	26	41	JIP 207101	5
83	27	62	JIP 207201	5
96	30	110	JIP 207301	5
110	37	160	JIP 207401	5
117	41	215	JIP 207421	5
127	42	240	JIP 207501	5
142	46	325	JIP 207521	5
154	50	420	JIP 207601	5
170	55	558	JIP 207701	5



EVAPORATION DISH, FLAT BOTTOM

with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
54	22	22	JIP 274101	5
70	30	62	JIP 274121	5
86	33	93	JIP 274201	5
98	40	154	JIP 274301	5
112	50	265	JIP 274401	5
126	53	330	JIP 274501	5
145	55	450	JIP 274601	5
151	63	620	JIP 274701	5
169	68	800	JIP 274801	1
195	74	1000	JIP 274901	1
254	105	2500	JIP 274991	1

