

VOLUMETRIC FLASKS

for sugar testing, with 2 marks and rim
borosilicate glass, calibrated "to contain" and marked at 20°C

Capacity	Cl. A	Pack of
50/55 ml	LMS 61319121	10
100/110 ml	LMS 61319122	10
200/220 ml	LMS 61319123	10



VOLUMETRIC FLASKS

for sugar testing acc. to Kohrausch with rim
borosilicate glass, calibrated "to contain" and marked at 20°C


Capacity	Cl. A	Pack of
100 ml	LMS 61313125	10
200 ml	LMS 61313126	10



SEPARATING FUNNELS

conical form (pear type) acc. to DIN ISO 4800, with plastic stopper,
graduation, borosilicate glass


Capacity	Neck size	PTFE	glass	Pack of
		stopcock key	stopcock key	
50 ml	NS 19/26	LMS 61165601	LMS 61165611	10
100 ml	NS 19/26	LMS 61165602	LMS 61165612	10
250 ml	NS 29/32	LMS 61165603	LMS 61165613	10
500 ml	NS 29/32	LMS 61165604	LMS 61165614	10
1000 ml	NS 29/32	LMS 61165605	LMS 61165615	2
2000 ml	NS 29/32	LMS 61165606	LMS 61165616	2



SEPARATING FUNNELS

acc. to Squibb acc. to DIN ISO 4800, with plastic stopper, graduation, borosilicate glass

Capacity	Neck size	PTFE	glass	Pack of
		stopcock key	stopcock key	
50 ml	NS 19/26	LMS 61165620	LMS 61165101	10
100 ml	NS 19/26	LMS 61165621	LMS 61165102	10
250 ml	NS 29/32	LMS 61165622	LMS 61165103	10
500 ml	NS 29/32	LMS 61165623	LMS 61165104	10
1000 ml	NS 29/32	LMS 61165624	LMS 61165105	2
2000 ml	NS 29/32	LMS 61165625	LMS 61165106	2



DROPPING FUNNELS

cylindrical form, acc. to DIN ISO 4800, with conical ground socket and cone,
plastic stopper, graduation, borosilicate glass

Capacity	NS size	PTFE	glass	Pack of
		stopcock key	stopcock key	
50 ml	NS 14/23	LMS 61164551	LMS 61164531	1
100 ml	NS 14/23	LMS 61164552	LMS 61164532	10
250 ml	NS 29/32	LMS 61164553	LMS 61164534	10
500 ml	NS 29/32	LMS 61164554	LMS 61164535	10
1000 ml	NS 29/32	LMS 61164555	LMS 61164536	10
2000 ml	NS 29/32	LMS 61164556	LMS 61164537	2

