


STAINING DISHES

acc. to Hellendahl, with wider upper section,
with cover, for 16 slides 76 x 26 mm, sodaglass


l x w x h	Code	Pack of
60 x 60 x 100 mm	HAV 23253111	10



STAINING DISHES

acc. to Schieferdecker,
with cover, for 10 slides 76 x 26 mm, sodaglass


l x w x h	Code	Pack of
90 x 70 x 40 mm	HAV 23253201	10



STAINING JAR

acc. to Coplin, with cover, for 10 slides 76 x 26 mm, sodaglass


w x h	Code	Pack of
40 / 100 mm	HAV 23253401	10



STAINING GLASS BOX

for staining tray, with cover, sodaglass


l x w x h	Code	Pack of
108 x 90 x 70 mm	HAV 23253423	10



STAINING TRAY

for 10 slides 76 x 26 mm, sodaglass

l x w x h	Code	Pack of
88 x 40 x 70 mm	HAV 23253422	10



STAINLESS STEEL HANDLE

for staining trays

Code	Pack of
HAV 23253427	10

